


EYETRONIC® - Restoring Vision in Glaucoma

41. Workshop Glaucoma, Carl Erb, MD, PhD
20th of January 2024




- Chronic stress leads to the termination of cell metabolism in the retinal ganglion cells and their nerve fibers.
- This leads to inactivity of the nerve cells so that their function – and Visual Field (VF) – are lost.
- Subsequently, nerve cells die causing deterioration of the optic nerve.



Optic Nerve Stimulation - Mechanism of Action

Preclinical evidence of the effects of electrical stimulation of the optic nerve


- Reactivation of nerve cells prevent further deterioration: Neuroprotection.
- Function is restored: Functional neurorestoration.
- Degenerated nerve fibers sprout and restore: Structural neurorestoration.



Optic Nerve Stimulation for patients with EYETRONIC®

External electrical stimulation with stimulation goggles leads to neuronal activity

- A biomarker confirms the evidence of signal transmission: the stimulation causes flashes of light, so called phosphenes.
- The cell metabolism of inactive nerve cells is reactivated and leads to the resumption of function.
- Neuronal signals are transmitted from the retina along the optic nerve to the visual cortex.



The diagram illustrates the EYETRONIC optic nerve stimulation process. It features four main components: 1) A pair of grey stimulation goggles with the 'eyetronic restoring vision' logo. 2) A close-up of the goggles' frame showing four red circular 'Embedded electrodes'. 3) A sagittal cross-section of a human head showing the brain, with a red lightning bolt indicating electrical stimulation and a red arrow labeled 'Gateway to the Brain' pointing from the eye to the brain. 4) A square image showing a bright blue starburst light labeled 'Biomarker' and 'Phosphenes'.

Painless, non-invasive stimulation of both eyes

CE certified Therapy

10 sessions for 40 - 70 minutes each

No serious adverse event in over 8 000 single sessions

EYETRONIC® - effectiveness in clinical routine use


More than 800 patients being treated with EYETRONIC® in clinical routine use

- The proof of therapy success is done via a standard perimetric examination of the visual field.
- Already three months after treatment, the visual field can improve.
- Impaired visual field areas can be restored in places.


Treatment with EYETRONIC®



Diagnosis with Standard Perimetry



Outcome: Improved Visual Field




Long-term data confirm EYETRONIC® efficacy


Key findings of 12-month PMCF study with 70 german glaucoma patients

- Perimetric examination after one year of 101 eyes in 70 glaucoma patients with progressive visual field loss despite IOP management.
- In 63.4% of treated eyes a halt in the progression of the visual field loss was observed.
- Improved Mean Defect (MD) and partial visual field restoration in 59%.


Eyes treated (n=101)



Comparison of Mean Defect (MD) (n=70)



Responder MD -2.2dB




Non Responder MD +0.7dB

Findings of 12-month follow-up of 7 Normal Tension Glaucoma (NTG) patients in 3 therapy centers


- Halt in visual field loss progression in 9 of 13 eyes (69.2%).
- Despite a small number of patients, significant improvement of visual field.

Halt of visual field loss in treated eyes (n=13)




■ Non Responder ■ Responder

MD of all patients @ 12 month



MD improved by -1.8 dB

MD change stronger in Responders



MD of responders improved by up to -8.2 dB

Peer-reviewed Journal Publication and presentation at the 10th World Glaucoma Congress (WGC)

12-month data in Bioelectronic Medicine, Erb et al., 2022

Normal Tension Glaucoma (NTG) patients, WGC 2023

Erb et al. *Bioelectronic Medicine* (2022) 8:6
<https://doi.org/10.1186/s42234-022-00089-9>

Bioelectronic Medicine

RESEARCH ARTICLE

Open Access

Electrical neurostimulation in glaucoma with progressive vision loss

Carl Erb¹, Sophie Eckert², Pia Gindorf¹, Martin Köhler³, Thomas Köhler³, Lukas Neuhann⁴, Thomas Neuhann⁴, Nadja Salzmann³, Stefanie Schmickler⁵ and Jens Ellrich^{6*}

Electrical optic nerve stimulation in normal tension glaucoma

P-303


J. ELLRICH¹, C. ERB², S. ECKERT³, M. KÖHLER⁴, T. KÖHLER⁴, N. SALZMANN⁴

email: glaucoma.ntg@gmail.com


¹ Medical Faculty, Friedrich-Alexander-University, Erlangen, Germany; ² Augenklinik Wittenbergplatz, Berlin, Germany; ³ Medizentrum, Neu-Ulm, Germany; ⁴ Gemeinschaftspraxis Salzmann/Köhler, Hannover, Germany

Right eye, NTG


PRE, 18.7 dB, 12 mmHg




2 w, 15 dB, 12 mmHg




3 m, 11.4 dB, 12 mmHg




6 m, 10.4 dB, 15 mmHg



9 m, 12.8 dB, 12 mmHg



12 m, 10.2 dB, 13 mmHg



Thank you very much for your attention – contact us!



E-Mail: therapie@eyetronic.com

E-Mail: therapy@eyetronic.com

www.eyetronic.com

Prevalence of glaucoma

- Crabb, Eye (Lond) 30: 304-13, 2016
- Esporcatte and Tavares, Arq Bras Oftalmol 79: 270-6, 2016
- Weinreb et al., JAMA 311:1901-11, 2014
- Quigley & Broman, Br J Ophthalmol 2006; 90:262–267, 2006

Glaucoma

- Cordeiro et al., BRAIN 2017: 140; 1757-1767
- De Moraes et al., Prog Ret Eye Res 56:107-47, 2017
- Wojcik-Gryciuk et al., Restor Neurol Neurosci 34: 107-23, 2015
- Garway-Heath et al, Lancet 2015; 385:1295-304
- Porciatti & Ventura, J Neuroophthalmol 32: 354-8, 2012
- Stanfield, Principles of Human Physiology, 4th Edition, Pearson, 2011
- Jones, Netter's Neurology, Elsevier, 2005
- Cordeiro et al., PNAS, vol 101, no 36, September 2004

Mode-of action of optic nerve stimulation

- Hanif et al., Exp Eye Res 149: 75-83, 2016
- Yin et al., Brain Res 1650: 10-20, 2016
- Fu et al., Graefes Arch Clin Exp Ophtalmol 253: 171-6, 2015
- Tagami et al., Jpn J Ophthalmol 53:257-66, 2009
- Miyake et al., Invest Ophthalmol Vis Sci 48:2356-61, 2007
- Morimoto et al., Invest ophthalmol Vis Sci 46: 2147-55, 2005

Mode-of action of optic nerve stimulation in Glaucoma

- A.H. Jassim et al., Ann Biomed Eng 2020, Sep 24, doi 10.1007/s10439

Clinical evidence of optic nerve stimulation

- Haberbosch et al., Front. Hum. Neurosci.13:43., 2019
- G Granata et al., Brain Stim 12 (2019), 800-802
- Erb, Ellrich, DOG, 2017
- Gall et al., PLOS ONE; 10.1371, 2016
- S Schmidt et al., Brian Stim 6 (2013), 87-93
- A Fedorov et al., Brain Stim 4 (2011), 189-201

Neuromodulation – Ophthalmology

- J Katz, Expert Review of Ophthalmology, 2016, Vol.11, No 5, 325-327
- Fu-Graefes, Arch Clin Exp Ophthalmol (2015) 253:171-176

Market

- Ophthalmology Equipment, MD 2334, marketsandmarkets.com, March 2018
- GlobalData „PharmaPoint: Glaucoma“, 2017
- ZionResearch: Stroke Diagnostics and Therapeutics Market, 2016
- Neuromodulation Market, marketsandmarkets.com, BT3571, July 2015
- Neuromodulation Market, modorintelligence.com, 2018
- Marketwatch, Neuromodulationmarket 2024, 2019